Technologies®  
**EMAG**  
www.emmi-dent.de

# EMMI-DENTAL HT / EMMI-DENT HT

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Technologies®  
**EMAG**  
www.emag-germany.de

EMAG AG  
DIN EN ISO 9001:2008

© EMAG AG  
2031111

**emmi-dental**  
 PROFESSIONAL

**emmi-dent 6**  
 PROFESSIONAL

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**emmi-dent 6**  
PROFESSIONAL

## It Functions so simply!

- 1** Connect the charging station with the socket.

- 2** Insert the handpiece from above into the charging station (charging time: ca. 24 hours)

- 3** Now insert the ultrasound attachment vertically from above into the handpiece and turn on the device.

- 4** Briefly wet the bristles of the ultrasound attachment under running water.

- 5** Apply only Ultrasound Toothpaste (the formation of millions of micro-bubbles only functions with this) to the bristles (approximately twice the size of a pea).

- 6** Lightly spread the toothpaste over the inside and outside of your teeth with the bristles of the ultrasound attachment.

- 7** Then hold the bristles motionless on the teeth, slightly from inside, outside and from above (chewing surface) – without applying pressure – for ca. 5 to 10 seconds. The bristles of the ultrasound attachment each cover 2-3 teeth. The actual retention time depends on the degree of impurity of the teeth.

- 8** After completion of dental cleaning (while the device is turned on), briefly rinse the bristles and handpiece under running water, turn the device off and briefly dry the handpiece (not the bristles). Ultrasound automatically disinfects the bristles.

- 9** You can simply plug the device with inserted ultrasound attachment into the charging station again after dental micro cleaning – or if you are for instance travelling, you can operate the device for ca. 20 days without recharging, if the batteries are fully charged (after ca. hours of charging time).

- 10** If you clean your teeth twice a day, please replace the ultrasound attachment after 3 months.

All products must be used only for the purposes and cases of application specified in the catalogues and technical descriptions, i.e. for cleaning teeth, and only in conjunction with nonoriginal devices and components recommended or approved by EMAG AG. Any other or uses above and beyond such uses are deemed improper use. EMAG AG does not accept any liability for any damage arising from improper use or guarantee the proper or perfect functioning of the device and its components where used improperly. Proper use also includes the observance of all of the instructions detailed in the operating instructions. This is furthermore indicated by the symbol  shown on the type plate. The perfect and safe operation of the products is subject to proper transport, proper storage, proper assembly and installation, as well as their careful operation and maintenance. Proper use also includes observance of the recurrent inspections detailed under the Accident Prevention Regulation A3.

The unit must be set up in such a way as to allow the charger's mains plug to be easily accessible at all times.

The non-conventional use of the product can be hazardous and cause damage to both life and limb, as well as to property. It is therefore prohibited to use the product for any uses other than its intended use. EMAG AG does not accept any liability for damage arising in connection with non-conventional use or guarantee the proper and perfect functioning of the product's components when used in non-conventional way.

- Only operate the ultrasonic dental micro cleaning and hygiene system at the line Voltage specified on the type plate (on the bottom of the charger)
- Do not clean the device with any flammable or aggressive chemicals (acids, etc.).
- In case of transport damages do not connect the charger to the mains.
- The device must only be operated in indoor spaces and at room temperature.

Remark: When switched on the green LED below the ON / OFF switch is only alight if the bristle attachment is plugged and if the ultrasound generator located in the bristle attachment works at the right frequency and correctly.

Battery must never be emptied completely ! Please pay attention that the hand piece is always switched off for re-charging because otherwise it will not be charged and may be damaged. Dentists recommend to apply cleaning of teeth twice a day for 3 to 5 minutes. If you clean your teeth twice a day, please exchange the ultrasonic head every 3 months. If you take the equipment with you when travelling the fully charged battery will be sufficient for approx. 20 days without recharging (if used twice a day for 2 x 4 minutes). Thereafter a re-charging for 24 hours is necessary. Otherwise the information of the original Operating Instructions are valid.

There is a guarantee of 24 months – from date of purchase – for this product. Defects which arise during this guarantee period in the form of material and/or manufacturing defects will be remedied at no charge, either by means of repair or through replacement delivery. Insofar as legally permissible, any other claims for compensation are excluded. The guarantee is inapplicable in the event of intervention by third parties or disassembly by third parties without our prior consent.

The guarantee is also inapplicable in the event of arbitrary damage or incorrect handling. The guarantee period is not extended with provision of a guarantee. We reserve the right to invoice the costs incurred for unjustified complaints such as installation or operating errors. For making a guarantee claim it is necessary to save the proof of purchase with date of purchase.

For fast handling of your guarantee claim cases, please send these – irrespective of specialised dealer – post-paid directly to EMAG AG. Unfortunately, unpaid shipments cannot be accepted due to logistical reasons.

Thank you for your cooperation and understanding.

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Date of Sale

Copyright EMAG AG. All rights reserved. Manufactured in FRG / EU. The contents of these operating instructions are subject to copyright and may neither be entirely or partially copied. Changes, abridgements, enhancements and additions, any publication or translation by third parties requires prior, written consent by EMAG AG.

EMAG AG  
DIN EN ISO 9001:2008

**RoHS** **CE**  
(§6Abs.3, ElektroG)

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© Copyright EMAG AG. Alle Rechte vorbehalten. Hergestellt in der BRD/EU. Der Inhalt dieser Bedienungsanleitung unterliegt dem Urheberrecht und sie darf weder ganz noch teilweise kopiert werden. Veränderungen, Kürzungen, Erweiterungen und Ergänzungen jede Veröffentlichung oder Übersetzung durch Dritte bedarf der vorherigen schriftlichen Einwilligung durch EMAG AG.